




ELEMISSION

LIBS TECHNOLOGY



APPLICATION NOTE



OPTIMIZING GOLD DEPOSIT EXPLORATION THROUGH RAPID
ELEMENTAL AND MINERALOGICAL CHARACTERIZATION USING ECORE





ELEMISION
LIBS TECHNOLOGY

Hf	Ta	W	Re	Os	Ir	Pt	Au	Hg	Tl
178.49	180.95	183.84	186.21	190.23	192.22	195.08	196.97	200.59	204.38
72	73	74	75	76	77	78	79	80	81
Zirconium	Tantalum	Tungsten	Rhenium	Osmium	Iridium	Platinum	Gold	Mercury	Thallium

Optimizing gold deposit exploration through rapid elemental and mineralogical characterization using ECORE

Importance of Gold Exploration

Gold exploration and production remains as a critical industry, driven by its economic importance (mining market valued at \$215.4 billion USD in 2024) and growing technological demand. Gold is consistently considered to be a trusted store of value, providing stability to financial markets, especially during periods of inflation or economic uncertainty. At the same time, technological advancements in the electronics, renewable energy, and aerospace industries have increased the demand for gold due to its superior conductivity, resistance to corrosion, and exceptional durability. The gold mining industry faces increasing challenges in sustaining production growth due to the diminishing availability of accessible gold deposits. The implementation of innovative technologies is essential to improve the discovery of new deposits and to optimize the recovery and production of existing deposits.

ECORE Technology



Figure 1: ELEMISION's ECORE Mobile laboratory, a LIBS-based core scanner mounted in a container for rapid onsite analysis of drill core.





Zr	Nb	Mo	Tc	Ru	Rh	Pd	Ag	Cu	Zn	Ga	Ge	As	Se	Br	Kr	Xe	Rn
91.22	92.91	95.94	[98]	101.07	102.91	106.42	107.87	112.41	114.82	118.71	121.76	127.60	132.91	148.91	171.87	186.21	222
Tungsten	Tungsten	Rhenium	Osmium	Iridium	Platinum	Gold	Mercury	Thallium	Lead	Bismuth	Po	Astatine	Polonium	Americium	Curium	Plutonium	Uranium
73	74	75	76	77	78	79	80	81	82	83	84	85	86	95	96	94	92
Hf	Ta	W	Re	Os	Ir	Pt	Au	Hg	Tl	Pb	Bi	Po	At	Am	Cm	Pu	U
178.49	180.95	183.84	186.21	190.23	192.22	195.08	196.97	200.59	204.38	207.2	208.98	209	210	227	238	238	238
Rutherfordium	Rhenium	Seaborgium	Rutherfordium	Rutherfordium	Rutherfordium	Rutherfordium	Rutherfordium	Rutherfordium	Rutherfordium	Rutherfordium	Rutherfordium	Rutherfordium	Rutherfordium	Rutherfordium	Rutherfordium	Rutherfordium	Rutherfordium
101	102	103	104	105	106	107	108	109	110	111	112	113	114	115	116	117	118

ECORE (Figure 1) utilizes Laser Induced Breakdown Spectroscopy (LIBS), interchangeably referred to as Laser Ablation Atomic Emission Spectroscopy (LA-AES). LIBS/LA-AES is a spectroscopic technique that utilizes laser ablation to excite elements in a material. Our ECORE technology employs this method while acquiring spectral data while scanning of drill core surfaces. With each laser pulse, a plasma is created, and as the plasma cools, the emitted light is captured by a spectrometer. This process generates a characteristic spectrum of each material ablated along the surface containing elements that are present in the material.

In this application note, we demonstrate how ECORE technology significantly accelerates gold deposit exploration by rapidly providing critical elemental and mineralogical information, greatly enhancing ore body knowledge.

APPLICATION: Distinguishing arsenopyrite from arsenic-bearing pyrite




Figure 2: A photograph, a mono-elemental arsenic (As) map, and a mineralogical map generated by ELEMISION's Smart Automated Mineralogy (SAM) algorithm of a section of drill core from an orogenic gold deposit. Arsenopyrite can be differentiated from As-bearing pyrite and non-As-bearing pyrite.



ELEMISSION
LIBS TECHNOLOGY

Zr	Nb	Mo	Tc	Ru	Rh	Pd	Ag	Cd	112.41	114
91.22	92.91	95.94	[98]	101.07	102.91	106.42	107.87	112.41	Mercury	Thall
72	73	74	75	76	77	78	79	80	8	
Hf	Ta	W	Re	Os	Ir	Pt	Au	Hg	T	
180.95	183.84	186.21	190.23	192.22	195.08	196.97	200.59			
Rhenium	Selenium	Bromine	Krypton							

For many gold deposit types, the presence of arsenic-bearing minerals is known to be associated with gold mineralization. Properly constraining the distribution of arsenic within a deposit can, therefore, provide invaluable insight into understanding controls on gold mineralization to facilitate decision-making and generate future drilling targets.

The orogenic gold deposit shown whose drill core is shown in Figure 2, has gold mineralization that is almost exclusively refractory and is associated with disseminated sulfide mineralization and related hydrothermal alteration. In general, gold particles are primarily trapped within fine-grained arsenopyrite or arsenic-rich pyrite crystals. Higher concentrations of arsenic are typically associated with higher gold concentrations.

The selectivity and sensitivity of ECORE technology allow users to distinguish between arsenopyrite, As-bearing pyrite, and non-As-bearing pyrite. Using a combination of mineralogical and elemental mapping, the distribution of arsenic throughout the core can be observed, and subsequently, mineralogical, textural, and chemical affinities between As-bearing pyrite and non-As-bearing pyrite can be established.

The textural and chemical characteristics of these minerals can be used to better understand the implications concerning the mechanisms and timing of gold deposition. Access to detailed mineralogy promotes easy and accurate deposit characterization and identification of alteration assemblages, allowing for informed decisions to be made for future exploration.

APPLICATION: Differentiating between carbonate types during exploration for quartz-carbonate vein-hosted gold deposits




Figure 3: A photograph, an elemental carbon (C) map, an elemental gold (Au) map, and a SAM image of a quartz-carbonate vein-hosted gold deposit drill core sample. Carbon-rich minerals and gold distribution can be clearly differentiated and visualized, along with their associated alteration minerals.





ELEMISION
LIBS TECHNOLOGY

Zr	Nb	Mo	Tc	Ru	Rh	Pd	Ag	Au	Hg	Pb
91.22	92.91	95.94	[98]	101.07	102.91	106.42	107.87	112.41	114	114
Hf	Ta	W	Re	Os	Ir	Pt	Au	Hg	Tl	Pb
178.49	180.95	183.84	186.21	190.23	192.22	195.08	196.97	200.59	204.38	207.2

Quartz-carbonate vein-hosted gold deposits typically exhibit mineralization closely associated with distinctive carbonate alteration assemblages, primarily comprising calcite and dolomite. Accurate identification and characterization of these carbonate minerals and related alteration phases is critical, as their presence and distribution often directly correlate with gold grades and the potential for economically viable mineralization.

In the example shown in Figure 3, gold mineralization is closely linked to quartz-carbonate veins and surrounding altered wall rocks, where gold particles occur both as free gold and in association with carbonate and sulfide minerals. Higher gold concentrations often correlate spatially with specific carbonate minerals, indicating a strong relationship between carbonate alteration assemblages and gold deposition.

ELEMISION's proprietary software for producing elemental and mineralogical (SAM) maps allows precise differentiation between carbon-rich minerals, gold occurrences, and other alteration minerals. Elemental maps, particularly those illustrating carbon and gold distributions, clearly shows spatial distribution and textural relationships of these minerals within the vein system.

Detailed mineralogical mapping and elemental analysis facilitate understanding of alteration zoning and its relationship to gold distribution. This enhanced characterization of carbonate alteration assemblages provides critical insights into the mineralizing processes, improves geological interpretations, and supports targeted exploration and resource delineation.





Zr	Nb	Mo	Tc	Ru	Rh	Pd	Ag	Au	Hg	Tl
91.22	92.91	95.94	[98]	101.07	102.91	106.42	107.87	112.41	200.59	204.38
Ta	W	Re	Os	Ir	Pt	Au	Hg	Tl	Pb	Bi
180.95	183.84	186.21	190.23	192.22	195.08	196.97	200.59	204.38	207.2	208.98
Seabornium	Roentgenium	Ubn	Uub	Uut	Uuq	Uup	Uuq	Uup	Uuh	Uuu
112	113	114	115	116	117	118	119	120	121	122

APPLICATION: Identifying base metal sulfide zonation as a vector for gold in VMS deposits





Figure 4: A photograph, elemental, maps, and a SAM image of a thin section offcut from a VMS deposit showing the association of base-metals to the gold present.



ELEMISSION

LIBS TECHNOLOGY



CONTACT INFORMATION

Website: www.elemission.ca

Email: info@elemission.ca

Phone: +1-574-998-3713

Address: 3410 Thimens Blvd, Saint Laurent, QC H4R 1V6

